

Anti-MAGP2 antibody ab47826

[1 References](#) [1 图像](#)

概述

产品名称	Anti-MAGP2抗体
描述	兔多克隆抗体to MAGP2
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Chimpanzee 
免疫原	Synthetic peptide corresponding to a region near the carboxy terminal end of MAGP2 protein (Human).
阳性对照	MAGP-2 transfected HEK293 lysate
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.01% Sodium azide Constituents: 0.42% Potassium phosphate, 0.87% Sodium chloride
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用 ab47826 于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度；实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		

应用说明

ELISA: 1/5000 - 1/20000.

IF (microscopy): Use at an assay dependent dilution.

WB: 1/250 - 1/1000. Predicted molecular weight: 19 kDa. .

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

靶标

功能

Component of the elastin-associated microfibrils.

序列相似性

Belongs to the MFAP family.


翻译后修饰

Forms intermolecular disulfide bonds either with other MAGP-2 molecules or with other components of the microfibrils.

细胞定位

Secreted > extracellular space > extracellular matrix.

图片



All lanes : Anti-MAGP2 antibody (ab47826) at 1/100 dilution

Lane 1 : Supernatants from non-transfected cells

Lane 2 : Secreted MAGP-2

Lane 3 : MAGP-2 transfected HEK293 lysate

Predicted band size: 19 kDa

Observed band size: 19 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors