

# Anti-MATN3/Matrilin-3 antibody ab106388

[2 References](#) [1 图像](#)

### 概述

---

产品名称	Anti-MATN3/Matrilin-3抗体
描述	兔多克隆抗体to MATN3/Matrilin-3
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Rat
免疫原	Synthetic peptide corresponding to Human MATN3/Matrilin-3 (internal sequence). 13 amino acid peptide from near the center of Human MATN3/Matrilin-3. Database link: <a href="#">AAI39908</a>
阳性对照	WB: Rat thymus tissue lysate.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 性能

---

形式	Liquid
存放说明	Shipped at 4°C. Store at 4°C (stable for up to 12 months).
存储溶液	pH: 7.2 Preservative: 0.02% Sodium azide Constituent: PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

### 应用

---

## The Abpromise guarantee

**Abpromise™** 承诺保证使用 ab106388 于以下的经测试应用


“应用说明”部分下显示的仅为推荐的起始稀释度；实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 53 kDa.

## 靶标

功能	Major component of the extracellular matrix of cartilage and may play a role in the formation of extracellular filamentous networks.
组织特异性	Expressed only in cartilaginous tissues, such as vertebrae, ribs and shoulders.
疾病相关	Multiple epiphyseal dysplasia 5 Spondyloepimetaphyseal dysplasia MATN3-related Osteoarthritis 2
序列相似性	Contains 4 EGF-like domains. Contains 1 VWFA domain.
细胞定位	Secreted.

## 图片



**Lane 1** : Anti-MATN3/Matrilin-3 antibody (ab106388) at 1 µg/ml

**Lane 2** : Anti-MATN3/Matrilin-3 antibody (ab106388) at 2 µg/ml

**All lanes** : Rat thymus tissue lysate

Lysates/proteins at 15 µg per lane.

**Predicted band size:** 53 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors