

Anti-MD1 antibody ab45424

2 References 2 图像

概述	
产品名称	Anti-MD1抗体
描述	兔多克隆抗体to MD1
宿主	Rabbit
经测试应用	适用于: IHC-P, WB
种属反应性	与反应: Human
免疫原	Synthetic 15 amino acid peptide from near the center of human MD1
阳性对照	Daudi cell lysate.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>
性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.2 Preservative: 0.02% Sodium azide Constituent: PBS
纯度	Immunogen affinity purified
纯化说明	This antibody is an immunoaffinity purified IgG
克隆	多克隆
同种型	IgG
应用	


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P		Use a concentration of 2.5 µg/ml.
WB		Use a concentration of 1 - 2 µg/ml. Detects a band of approximately 27 kDa (predicted molecular weight: 18 kDa).

靶标

功能	May cooperate with CD180 and TLR4 to mediate the innate immune response to bacterial lipopolysaccharide (LPS) and cytokine production. Important for efficient CD180 cell surface expression.
组织特异性	Highly expressed in B-cells, monocytes and tonsil.
细胞定位	Secreted > extracellular space. Associated with CD180 at the cell surface.


图片



Anti-MD1 antibody (ab45424) at 2 µg/ml + Daudi cell lysate at 15 µg

Predicted band size: 18 kDa

Observed band size: 27 kDa



ab45424 at 2.5µg/ml staining MD1 in human spleen tissue section by Immunohistochemistry (Formalin/ PFA fixed paraffin-embedded tissue sections).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors