

Anti-Mephalan antibody [MP5/73] ab97892

3 References

概述

产品名称	Anti-Mephalan抗体[MP5/73]
描述	大鼠单克隆抗体[MP5/73] to Mephalan
宿主	Rat
特异性	This antibody was developed to improve detection of melphalan-DNA adducts through its ability to recognise the stabilised adduct. This antibody can be used for quantification of melphalan adducts in purified DNA by ELISA and for immunofluorescent detection of melphalan-DNA adducts in single cells.
经测试应用	适用于: ELISA
种属反应性	与反应: Species independent
免疫原	Full length native protein (purified) corresponding to Mephalan.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.02% Sodium azide Constituent: 99.98% PBS
纯度	Protein A/G purified
克隆	单克隆
克隆编号	MP5/73
骨髓瘤	Y3/Ag1.2.3
同种型	IgG2b

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab97892于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		Use at an assay dependent dilution.

靶标

相关性 Mephalan hydrochloride is a chemotherapy drug belonging to the class of nitrogen mustard alkylating agents.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors