

Anti-Met (c-Met) antibody [EP1454Y] - BSA and Azide free ab157370

敲除验证
重组
RabMAb

7 图像

概述

产品名称	Anti-Met (c-Met)抗体[EP1454Y] - BSA and Azide free
描述	兔单克隆抗体[EP1454Y] to Met (c-Met) - BSA and Azide free
宿主	Rabbit
经测试应用	<p>适用于: IHC-P, WB, Indirect ELISA</p> <p>不适用于: Flow Cyt or ICC/IF</p>
种属反应性	与反应: Mouse, Rat, Human
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
阳性对照	WB: Wild-type HAP1 cell lysate. HepG2, HEK-293 and HeLa cell lysate. Mouse and rat thymus tissue lysate. Mouse lung tissue lysate. Hela and A459 whole cell lysate. IHC-P: Human bladder carcinoma and clear cell kidney carcinoma tissue.
常规说明	<p>ab157370 is the carrier-free version of ab51067.</p> <p>Our carrier-free antibodies are typically supplied in a PBS-only formulation, purified and free of BSA, sodium azide and glycerol. The carrier-free buffer and high concentration allow for increased conjugation efficiency.</p> <p>This conjugation-ready format is designed for use with fluorochromes, metal isotopes, oligonucleotides, and enzymes, which makes them ideal for antibody labelling, functional and cell-based assays, flow-based assays (e.g. mass cytometry) and Multiplex Imaging applications.</p> <p>Use our conjugation kits for antibody conjugates that are ready-to-use in as little as 20 minutes with <1 minute hands-on-time and 100% antibody recovery: available for fluorescent dyes, HRP, biotin and gold.</p> <p>This product is compatible with the Maxpar[®] Antibody Labeling Kit from Fluidigm, without the need for antibody preparation. Maxpar[®] is a trademark of Fluidigm Canada Inc.</p> <p>Our RabMAb[®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents.</p>

性能

形式	Liquid
----	--------

存放说明	Shipped at 4°C. Store at +4°C. Do Not Freeze.
存储溶液	Constituent: PBS
无载体	是
纯度	Protein A purified
纯化说明	Protein-A purification via MabSelect SuRe
克隆	单克隆
克隆编号	EP1454Y
同种型	IgG

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab157370于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P		Use at an assay dependent concentration. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.
WB		Use at an assay dependent concentration. Detects a band of approximately 160 kDa (predicted molecular weight: 156 kDa).
Indirect ELISA		Use at an assay dependent concentration.

应用说明 Is unsuitable for Flow Cyt or ICC/IF.

靶标

功能	Receptor for hepatocyte growth factor and scatter factor. Has a tyrosine-protein kinase activity. Functions in cell proliferation, scattering, morphogenesis and survival.
疾病相关	<p>Note=Activation of MET after rearrangement with the TPR gene produces an oncogenic protein.</p> <p>Note=Defects in MET may be associated with gastric cancer.</p> <p>Defects in MET are a cause of hepatocellular carcinoma (HCC) [MIM:114550].</p> <p>Defects in MET are a cause of renal cell carcinoma papillary (RCCP) [MIM:605074]. It is a subtype of renal cell carcinoma tending to show a tubulo-papillary architecture formed by numerous, irregular, finger-like projections of connective tissue. Renal cell carcinoma is a heterogeneous group of sporadic or hereditary carcinoma derived from cells of the proximal renal tubular epithelium. It is subclassified into common renal cell carcinoma (clear cell, non-papillary carcinoma), papillary renal cell carcinoma, chromophobe renal cell carcinoma, collecting duct carcinoma with medullary carcinoma of the kidney, and unclassified renal cell carcinoma.</p> <p>Note=A common allele in the promoter region of the MET shows genetic association with susceptibility to autism in some families. Functional assays indicate a decrease in MET promoter activity and altered binding of specific transcription factor complexes.</p> <p>Note=MET activating mutations may be involved in the development of a highly malignant, metastatic syndrome known as cancer of unknown primary origin (CUP) or primary occult malignancy. Systemic neoplastic spread is generally a late event in cancer progression. However,</p>

in some instances, distant dissemination arises at a very early stage, so that metastases reach clinical relevance before primary lesions. Sometimes, the primary lesions cannot be identified in spite of the progresses in the diagnosis of malignancies.

序列相似性

Belongs to the protein kinase superfamily. Tyr protein kinase family.
Contains 3 IPT/TIG domains.
Contains 1 protein kinase domain.
Contains 1 Sema domain.

结构域

The kinase domain is involved in SPSB1 binding.


翻译后修饰

Dephosphorylated by PTPRJ at Tyr-1349 and Tyr-1365.

细胞定位

Membrane.

图片



All lanes : Anti-Met (c-Met) antibody [EP1454Y] - N-terminal
([ab51067](#)) at 1/1000 dilution

Lane 1 : Wild-type A549 cell lysate

Lane 2 : MET knockout A549 cell lysate

Lane 3 : Wild-type HeLa [ab255929](#) cell lysate

Lane 4 : MET (Met (c-Met)) knockout HeLa [ab256991](#) cell lysate

Lane 5 : HT-29 cell lysate

Lane 6 : T-47D cell lysate

Lysates/proteins at 20 µg per lane.

Developed using the ECL technique.

Performed under reducing conditions.


Predicted band size: 156 kDa

Observed band size: 40-50 kDa

This data was developed using the same antibody clone in a different buffer formulation ([ab51067](#)).

Western blot: Anti-MET antibody [EP1454Y] ([ab51067](#)) staining at 1/1000 dilution, shown in green; Mouse anti-GAPDH antibody [6C5] ([ab8245](#)) loading control staining at 1/20000 dilution, shown in magenta. In Western blot, [ab51067](#) was shown to bind specifically to MET. A band was observed at 40-50 kDa in wild-type A549 cell lysates with no signal observed at this size in MET knockout cell line. To generate this image, wild-type and MET knockout A549 cell lysates were analysed. First, samples were run on an SDS-PAGE

gel then transferred onto a nitrocellulose membrane. Membranes were blocked in 3% milk in TBS-0.1% Tween® 20 (TBS-T) before incubation with primary antibodies overnight at 4°C. Blots were washed four times in TBS-T, incubated with secondary antibodies for 1 h at room temperature, washed again four times then imaged. Secondary antibodies used were Goat anti-Rabbit IgG H&L 800CW and Goat anti-Mouse IgG H&L 680RD at 1/20000 dilution.



All lanes : Anti-Met (c-Met) antibody [EP1454Y] - N-terminal ([ab51067](#)) at 1/1000 dilution

All lanes :

Lysates/proteins at 20 µg per lane.


Performed under reducing conditions.

Predicted band size: 156 kDa

Observed band size: 50 kDa

This data was developed using the same antibody clone in a different buffer formulation ([ab51067](#)).

False colour image of Western blot: Anti-Met (c-Met) antibody [EP1454Y] - N-terminal staining at 1/1000 dilution, shown in green; Rabbit Anti-GAPDH antibody [EPR16891] ([ab181602](#)) loading control staining at 1/20000 dilution, shown in red. In Western blot, [ab51067](#) was shown to bind specifically to the alpha chain of c-Met. A band was observed at 50 kDa in wild-type HeLa cell lysates with no signal observed at this size in MET knockout cell line [ab265961](#) (knockout cell lysate [ab256991](#)). To generate this image, wild-type and MET knockout HeLa cell lysates were analysed. First, samples were run on an SDS-PAGE gel then transferred onto a nitrocellulose membrane. Membranes were blocked in 3 % milk in TBS-0.1 % Tween® 20 (TBS-T) before incubation with primary antibodies overnight at 4°C. Blots were washed four times in TBS-T, incubated with secondary antibodies for 1 h at room temperature, washed again four times then imaged. Secondary antibodies used were HRP conjugated Goat anti-Rabbit (H+L) and Goat anti-Mouse IgG H&L (IRDye® 680RD) preabsorbed ([ab216776](#)) at 1/20000 dilution.




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Met (c-Met) antibody [EP1454Y] - BSA and Azide free (ab157370)

Immunohistochemical staining of paraffin embedded human clear cell kidney carcinoma with purified [ab51067](#) at a working dilution of 1/100. The secondary antibody used is HRP goat anti-rabbit IgG H&L ([ab97051](#)) at 1/500. The sample is counterstained with hematoxylin. Antigen retrieval was performed using Tris-EDTA buffer, pH 9.0.

PBS was used instead of the primary antibody as the negative control (inset).

This data was developed using the same antibody clone in a different buffer formulation containing PBS, BSA, glycerol, and sodium azide ([ab51067](#)).




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Met (c-Met) antibody [EP1454Y] - BSA and Azide free (ab157370)

This image is courtesy of an Abreview submitted by David Ivancic.

[ab51067](#) staining Met (c-Met) in human breast tissue sections by Immunohistochemistry (IHC-P - paraformaldehyde-fixed, paraffin-embedded sections).

Tissue was fixed with formaldehyde, permeabilized with 0.05% Tween-20 and blocked for 30 minutes at 22°C; antigen retrieval was by heat mediation in antigen retrieval buffer (100X citrate buffer pH 6.0) ([ab94674](#)). Samples were incubated with the primary antibody (1/100) for 14 hours at 4°C. An Alexa Fluor® 488-conjugated goat anti-rabbit IgG polyclonal (1/300) was used as the secondary antibody.

This data was developed using the same antibody clone in a different buffer formulation containing PBS, BSA, glycerol, and sodium azide ([ab51067](#)).




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Met (c-Met) antibody [EP1454Y] - BSA and Azide free (ab157370)

Immunohistochemical staining of paraffin embedded human bladder carcinoma with purified **ab51067** at a working dilution of 1/100. The secondary antibody used is HRP goat anti-rabbit IgG H&L (**ab97051**) at 1/500. The sample is counterstained with hematoxylin. Antigen retrieval was performed using Tris-EDTA buffer, pH 9.0.

PBS was used instead of the primary antibody as the negative control (inset).

This data was developed using the same antibody clone in a different buffer formulation containing PBS, BSA, glycerol, and sodium azide (**ab51067**).

Indirect ELISA antibody dose-response curve antigen at 1000 ng/ml




Indirect ELISA - Anti-Met (c-Met) antibody [EP1454Y] - BSA and Azide free (ab157370)

This data was developed using **ab51067**, the same antibody clone in a different buffer formulation.

ELISA analysis of Human Met recombinant protein at 1000 ng/mL with **ab51067**. An Alkaline Phosphatase-conjugated AffiniPure Goat Anti-Rabbit IgG (H+L) at 1/2500 dilution was used as the secondary antibody.

Why choose a recombinant antibody?



Anti-Met (c-Met) antibody [EP1454Y] - BSA and Azide free (ab157370)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors