

Anti-Metapneumovirus Phosphoprotein antibody [hMPV 33] ab94801

3 References

概述

产品名称	Anti-Metapneumovirus Phosphoprotein抗体[hMPV 33]
描述	小鼠单克隆抗体[hMPV 33] to Metapneumovirus Phosphoprotein
宿主	Mouse
特异性	Identification and location of expression of human MPV phosphoprotein of both types A and B.
经测试应用	适用于: ELISA, ICC/IF
种属反应性	与反应: Metapneumovirus
免疫原	Tissue, cells or virus corresponding to Metapneumovirus Phosphoprotein. Human hMPV strain NCL-174
阳性对照	Human MPV infected LLC-MK2 cells
常规说明	<p>Fusion partner: PS-NS/1-Ag4 Growth medium RPMI with 25mM Hepes, 2mM additional glutamine and 10% heat inactivated foetal calf serum</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	<p>Preservative: 0.02% Sodium azide</p> <p>Constituent: 99.98% PBS</p>
纯度	Protein A purified
克隆	单克隆
克隆编号	hMPV 33

骨髓瘤 P3-NS1/1-Ag4-1
同种型 IgG2a

应用

The Abpromise guarantee Abpromise™承诺保证使用ab94801于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		Use at an assay dependent concentration.
ICC/IF		Use at an assay dependent concentration. Acetone fixed

靶标

相关性 Human metapneumovirus (hMPV) is a member of the Paramyxoviridae family belonging to the Metapneumovirus genus.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors