

Anti-Mitofilin antibody ab93323

[3 References](#) [3 图像](#)

概述

产品名称	Anti-Mitofilin抗体
描述	兔多克隆抗体to Mitofilin
宿主	Rabbit
经测试应用	适用于: ICC/IF
种属反应性	与反应: Human
免疫原	Recombinant full length protein (Human)
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.05% Sodium azide Constituent: Whole serum
纯度	Whole antiserum
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab93323于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
ICC/IF		1/500.

靶标

细胞定位


Mitochondrion inner membrane.

图片




Immunofluorescent analysis of Mitofilin using anti-Mitofilin polyclonal antibody (ab93323) shows staining in p19 Cells.

Immunocytochemistry/ Immunofluorescence - Anti-Mitofilin antibody (ab93323)



Immunofluorescent analysis of Mitofilin using anti-Mitofilin polyclonal antibody (ab93323) shows staining in HMVEC Cells.

Immunocytochemistry/ Immunofluorescence - Anti-Mitofilin antibody (ab93323)



Immunofluorescent analysis of Mitofilin using anti-Mitofilin polyclonal antibody (ab93323) shows staining in A549 Cells.

Immunocytochemistry/ Immunofluorescence - Anti-Mitofilin antibody (ab93323)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors