

Anti-mono methyl Arginine antibody [16B11] ab414

14 References

概述

产品名称	Anti-mono methyl Arginine抗体[16B11]
描述	小鼠单克隆抗体[16B11] to mono methyl Arginine
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Species independent
常规说明	

This antibody will be of central importance in analysing the methylation status of chromatin and transcription factors. Since it recognises monomethyl arginine only it can be used in parallel with [ab413](#) (detects dimethyl arginine only) and [ab412](#) (detects both monomethyl and dimethyl arginine) to monitor the exact modification status.

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Constituents: 1.21% Tris, 0.75% Glycine, 2% Sucrose
纯度	Protein A purified
Primary antibody说明	This antibody will be of central importance in analysing the methylation status of chromatin and transcription factors. Since it recognises monomethyl arginine only it can be used in parallel with ab413 (detects dimethyl arginine only) and ab412 (detects both monomethyl and dimethyl arginine) to monitor the exact modification status.
克隆	单克隆
克隆编号	16B11

骨髓瘤	unknown
同种型	IgG2b
轻链类型	kappa

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab414于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度；实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/50 - 1/500.

靶标

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors