

Anti-Mouse IgE antibody ab19967

2 References

概述

产品名称	抗小鼠IgE抗体
描述	山羊多克隆抗体to小鼠IgE
宿主	Goat
特异性	Reacts only with IgE by IEP, gel diffusion and ELISA.
经测试应用	适用于: Immunoelectrophoresis, ELISA
种属反应性	与反应: Mouse
免疫原	Full length protein corresponding to Mouse Mouse IgE. Native protein
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at +4°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.1% Sodium azide Constituent: PBS
纯度	Immunogen affinity purified
纯化说明	Immunoaffinity purified off an immunoabsorbant containing a mouse IgE protein different from that used for immunization.
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™承诺保证使用ab19967于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Immunoelectrophoresis		Use at an assay dependent concentration.
ELISA		Use at an assay dependent concentration.

靶标

相关性 IgE mediates the immune responses in the pathogenesis of allergic asthmas. This has been shown in the successful use of recombinant 'humanised' monoclonal antibodies against IgE in the treatment of moderate-to-severe allergic asthma where the antibody forms complexes with free IgE and blocks its interaction with mast cells and basophils.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors