

Mouse IgG1, kappa monoclonal [MG1-45] - Isotype control ab18448

9 References

概述

产品名称	小鼠IgG1, kappa单克隆抗体[MG1-45]-同型对照
经测试应用	适用于: IHC-P, IHC-Fr, Flow Cyt, WB, IP, ELISA
免疫原	Chemical/ Small Molecule (Trinitrophenol) conjugated to KLH.
常规说明	<p>The ab18448, MG1-45 immunoglobulin reacts with KLH and was chosen as an isotype control after screening on a variety of resting, activated, live, and fixed mouse and human tissues.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.20 Preservative: 0.09% Sodium azide Constituent: PBS
纯度	Protein G purified
纯化说明	The immunoglobulin was purified by affinity chromatography.
Isotype control说明	The ab18448, MG1-45 immunoglobulin reacts with KLH and was chosen as an isotype control after screening on a variety of resting, activated, live, and fixed mouse and human tissues.
克隆	单克隆
克隆编号	MG1-45
同种型	IgG1

应用

The Abpromise guarantee

Abpromise™承诺保证使用ab18448于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.
Flow Cyt		Use at an assay dependent concentration. Suitable for surface or intracellular staining.
WB		Use at an assay dependent concentration.
IP		Use at an assay dependent concentration.
ELISA		Use at an assay dependent concentration. PubMed: 18006877

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors