

### Mouse IgG3, monoclonal [PPV-07] - Isotype Control ab91537

#### 概述

产品名称	小鼠IgG3,单克隆抗体[PPV-07] -同型对照
特异性	ab91537 reacts with an undefined epitope on a plant pathogen. No reactivity with cell surface antigens on Human, Porcine, Mouse or Rat. It can be used as a negative control for immunofluorescence staining (e.g. Flow Cytometry) when using Mouse monoclonal antibodies to Human cell surface antigens of isotype IgG1. ab91537 enables an estimation of non specific binding of Mouse monoclonal antibodies to cell surface components in peripheral blood and tissue. It is suitable for whole blood, Ficoll separated preparations, frozen and paraffin embedded sections.
经测试应用	适用于: Flow Cyt
免疫原	not applicable
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	<p>pH: 7.20</p> <p>Preservative: 0.09% Sodium azide</p> <p>Constituent: 1% BSA</p>
纯度	Immunogen affinity purified
纯化说明	Purified (pass through 0.22 µm filter).
克隆	单克隆
克隆编号	PPV-07
同种型	IgG3

## 应用

---

**The Abpromise guarantee**      **Abpromise™**承诺保证使用ab91537于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		

### 应用说明

Flow Cyt: Use 10 µl for 10<sup>6</sup> cells.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors