

Anti-MPP8 antibody ab122611

3 图像

概述	
产品名称	Anti-MPP8抗体
描述	兔多克隆抗体to MPP8
宿主	Rabbit
经测试应用	适用于: IHC-P, WB, ICC/IF
种属反应性	与反应: Human
免疫原	antigen sequence, corresponding to amino acids 421-528 of Human MPP8.
阳性对照	Human Colon tissue; RT 4, U 251 MG and Human Tonsil lysates; U 251 MG cells.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>
性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at -20°C or -80°C. Avoid freeze / thaw cycle. Store undiluted.
存储溶液	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol (glycerin, glycerine)
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG
应用	

The Abpromise guarantee **Abpromise™**承诺保证使用ab122611于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P		1/1000 - 1/2500. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.
WB		1/250 - 1/500.
ICC/IF		Use a concentration of 0.25 - 2 µg/ml.

靶标

序列相似性	Contains 4 ANK repeats. Contains 1 chromo domain.
翻译后修饰	Phosphorylated in M (mitotic) phase.
细胞定位	Nucleus.


图片



All lanes : Anti-MPP8 antibody (ab122611) at 1/250 dilution


- Lane 1 :** RT 4 lysate
- Lane 2 :** U 251 MG lysate
- Lane 3 :** Human Plasma lysate
- Lane 4 :** Human Liver lysate
- Lane 5 :** Human Tonsil lysate

Developed using the ECL technique.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-MPP8 antibody (ab122611)

ab122611, at 1/1000 dilution, staining MPP8 in Paraffin-embedded Human Colon tissue by Immunohistochemistry.



Immunocytochemistry/ Immunofluorescence - Anti-MPP8 antibody (ab122611)

ab122611, at 4µg/ml, staining MPP8 in PFA/Triton X-100 fixed/permealised U 251 MG cells by Immunofluorescence.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors