

Anti-MTMR3 antibody ab122213

2 图像

| | |
|-------|---|
| 概述 | |
| 产品名称 | Anti-MTMR3抗体 |
| 描述 | 兔多克隆抗体to MTMR3 |
| 宿主 | Rabbit |
| 经测试应用 | 适用于: IHC-P, ICC/IF |
| 种属反应性 | 与反应: Human |
| 免疫原 | antigen sequence, corresponding to amino acids 636-740 of Human MTMR3. |
| 阳性对照 | Human testis tissue; A431 cells. |
| 常规说明 | <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p> |
| 性能 | |
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |
| 存储溶液 | pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 40% Glycerol (glycerin, glycerine), 59% PBS |
| 纯度 | Immunogen affinity purified |
| 克隆 | 多克隆 |
| 同种型 | IgG |
| 应用 | |

The Abpromise guarantee

Abpromise™承诺保证使用ab122213于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

| 应用 | Ab评论 | 说明 |
|--------|------|--|
| IHC-P | | 1/1000. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol. |
| ICC/IF | | Use a concentration of 1 - 4 µg/ml. Fixation/Permeabilization: PFA/Triton X-100 |

靶标


| | |
|-------|--|
| 功能 | Phosphatase that acts on lipids with a phosphoinositol headgroup. Has phosphatase activity towards phosphatidylinositol-3-phosphate and phosphatidylinositol-3,5-bisphosphate. May also dephosphorylate proteins phosphorylated on Ser, Thr, and Tyr residues. |
| 序列相似性 | Belongs to the protein-tyrosine phosphatase family. Non-receptor class myotubularin subfamily. Contains 1 FYVE-type zinc finger. Contains 1 myotubularin phosphatase domain. |
| 细胞定位 | Cytoplasm. Membrane. |

图片



ab122213 at 1/1000 dilution staining MTMR3 in Paraffin-embedded Human Testis tissue by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-MTMR3 antibody (ab122213)



ab122213 at 4 µg/ml staining MTMR3 (green) in PFA-fixed/Triton-permeabilized Human A431 cells by Immunofluorescence.

Immunocytochemistry/ Immunofluorescence - Anti-MTMR3 antibody (ab122213)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors