



Anti-Mumps nucleoprotein antibody [8H4] ab9876

7 References **1 图像**

概述

产品名称	Anti-Mumps nucleoprotein抗体[8H4]
描述	小鼠单克隆抗体[8H4] to Mumps nucleoprotein
宿主	Mouse
特异性	Recognises recombinant yeast-derived mumps NP. Recognises native NP in mumps virus infected Vero cells. This antibody has been tested on 3 different mumps virus strains isolated in the UK: 55799-98, SJ801-95 and Gloucester.
经测试应用	适用于: ICC/IF
种属反应性	与反应: Mumps virus
免疫原	Full length protein corresponding to Mumps nucleoprotein aa 1 to the C-terminus. Database link: AAG37826.1
常规说明	<div>  Run BLAST with  Run BLAST with </div> <p>This product was changed from ascites to tissue culture supernatant on 28/11/2017. Please note that the dilutions may need to be adjusted accordingly.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 6.80 Preservative: 0.1% Sodium azide Constituent: 15% Fetal Bovine Serum (FBS)
纯度	Tissue culture supernatant
克隆	单克隆

克隆编号	8H4
骨髓瘤	NS0/1
同种型	IgG2b
轻链类型	unknown

应用


The Abpromise guarantee **Abpromise™**承诺保证使用ab9876于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		1/1 - 1/20.

靶标

图片



Immunofluorescence staining of mumps virus infected Vero cells with monoclonal antibody 8H4.

Immunofluorescence staining of mumps virus infected Vero cells with monoclonal antibody 8H4.

Immunocytochemistry/ Immunofluorescence - Anti-Mumps nucleoprotein antibody [8H4] (ab9876)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors