

# Anti-Mycoplasma bovis antibody [MY11-320.1] ab62385

## 1 References

### 概述

产品名称	Anti-Mycoplasma bovis抗体[MY11-320.1]
描述	小鼠单克隆抗体[MY11-320.1] to Mycoplasma bovis
宿主	Mouse
经测试应用	适用于: ELISA
种属反应性	与反应: Mycoplasma bovis
免疫原	Tissue, cells or virus corresponding to Mycoplasma bovis. Mycoplasma bovis cells.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.40 Constituent: PBS
纯度	Protein G purified
克隆	单克隆
克隆编号	MY11-320.1
同种型	IgG1

### 应用

**The Abpromise guarantee**      **Abpromise™** 承诺保证使用ab62385于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		Use at an assay dependent concentration.

## 靶标

**相关性** Since it was first recognized in the early 1960's, mycoplasma mastitis has been reported in many dairy farms throughout the United States. Mycoplasma are opportunistic invaders and can cause joint infections, pneumonia, metritis, and mastitis. Several species of Mycoplasma have been found in dairy cattle, but by far the most prevalent is Mycoplasma bovis, which is responsible for more than half of all reported isolates. This specie also causes the most severe infections of any Mycoplasma organism known to cause disease in cattle.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors