

Anti-Myofibroblasts antibody [PR 2D3] ab20255

★★★★★ [1 Abreviews](#) [2 References](#)

概述

产品名称	Anti-Myofibroblasts抗体[PR 2D3]
描述	小鼠单克隆抗体[PR 2D3] to Myofibroblasts
宿主	Mouse
经测试应用	适用于: IHC-Fr
种属反应性	与反应: Human
免疫原	Tissue, cells or virus corresponding to Human Myofibroblasts. normal colorectal mucosa (Human).
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.02% Sodium azide Constituent: 99.98% PBS
纯度	Protein A/G purified
克隆	单克隆
克隆编号	PR 2D3
骨髓瘤	NS0
同种型	IgG1
轻链类型	unknown

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab20255于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-Fr		

应用说明

IHC-Fr, Acetone Fixed: Use at an assay dependent dilution.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

靶标

相关性

The antibody reacts with a cell membrane component of cells in the pericrypt sheath, with smooth muscle and with myofibroblasts. The gold standard for identification of myofibroblasts.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors