


Anti-Myogenin antibody ab103924

★★★★★ [1 Abreviews](#) [2 References](#) [1 图像](#)

概述

产品名称	Anti-Myogenin抗体
描述	山羊多克隆抗体to Myogenin
宿主	Goat
经测试应用	适用于: WB
种属反应性	与反应: Mouse 预测可用于: Rat, Sheep, Rabbit, Goat, Cow, Dog, Human, Pig 
免疫原	Synthetic peptide: C-TDAHNLHSLTSIVD , corresponding to C terminal amino acids 193-206 of Human Myogenin (NP_002470.2). Run BLAST with Run BLAST with
阳性对照	Recombinant Human Myogenin protein (ab114244) can be used as a positive control in WB. Mouse Skeletal Muscle lysate
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	pH: 7.30 Preservative: 0.02% Sodium azide Constituents: Tris buffered saline, 0.5% BSA
纯度	Immunogen affinity purified
纯化说明	ab103924 was purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

克隆 多克隆
同种型 IgG

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab103924于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		Use a concentration of 2 - 4 µg/ml. Detects a band of approximately 38 kDa (predicted molecular weight: 25 kDa). A 1 hour primary incubation is recommended for this product.


靶标

功能 Involved in muscle differentiation (myogenic factor). Induces fibroblasts to differentiate into myoblasts. Probable sequence specific DNA-binding protein.

序列相似性 Contains 1 basic helix-loop-helix (bHLH) domain.

细胞定位 Nucleus.

图片



Anti-Myogenin antibody (ab103924) at 2 µg/ml + Mouse Skeletal Muscle lysate in RIPA buffer at 35 µg

Developed using the ECL technique.

Predicted band size: 25 kDa

Observed band size: ~38 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors