

Anti-N-epsilon acetyl Lysine antibody [11A1] ab409

9 References 1 图像


概述	
产品名称	Anti-N-epsilon acetyl Lysine抗体[11A1]
描述	小鼠单克隆抗体[11A1] to N-epsilon acetyl Lysine
宿主	Mouse
经测试应用	适用于: WB, ELISA, IHC-P
种属反应性	与反应: Species independent
常规说明	<p>The ability of this antibody to fix complement is unknown.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>
性能	
形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Constituents: 0.75% Glycine, 1.21% Tris, 2% Sucrose
纯度	Protein A purified
克隆	单克隆
克隆编号	11A1
骨髓瘤	unknown
同种型	IgG1
轻链类型	kappa
应用	

The Abpromise guarantee **Abpromise™**承诺保证使用ab409于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/500 - 1/5000.
ELISA		1/5000 - 1/50000.
IHC-P		1/500 - 1/5000.

图片



Western Blot analysis showing ab409 staining N-epsilon acetyl Lysine on immunoprecipitated ER-resident proteins from control (non-transfected) H4 cells.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors