

# Anti-NAS6 antibody ab91447

**2 References**   [1 图像](#)

### 概述

产品名称	Anti-NAS6抗体
描述	兔多克隆抗体to NAS6
宿主	Rabbit
经测试应用	适用于: WB, IP
种属反应性	与反应: <i>Saccharomyces cerevisiae</i>
免疫原	Recombinant full length <i>Saccharomyces cerevisiae</i> NAS6 protein, expressed in <i>E. coli</i> .
阳性对照	Crude extract of <i>S. cerevisiae</i>
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	<p>pH: 6</p> <p>Preservative: 0.09% Sodium azide</p> <p>Constituents: PBS, 0.1% BSA, 50% Glycerol (glycerin, glycerine)</p>
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

### 应用

**The Abpromise guarantee**   **Abpromise™**承诺保证使用ab91447于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
WB		1/5000. Predicted molecular weight: 26 kDa.
IP		Use at an assay dependent dilution.

## 靶标

**相关性** NAS6 (p28) is a proteasome-interacting protein which acts as a chaperone to assemble the base sub-complex of the 19S RP which is involved in the ATP-dependent degradation of ubiquitinated proteins. NAS6 is composed of 226 amino acid residues with calculated molecular masses of 24.4kD.

**细胞定位** Cytosol; Nucleus; Proteasome complex

## 图片



Anti-NAS6 antibody (ab91447) at 1/5000 dilution + Crude extract of *S. cerevisiae*

**Predicted band size: 26 kDa**

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors