

### Anti-NCAM1 antibody [123C3] ab9272

★★★★ 4 Abreviews 23 References 6 图像

#### 概述

产品名称	Anti-NCAM1抗体[123C3]
描述	小鼠单克隆抗体[123C3] to NCAM1
宿主	Mouse
特异性	This antibody stains neural cell adhesion molecule (NCAM) which is present in small cell lung cancer and lung carcinoids. It also reacts with a number of non-small cell lung carcinomas and neuroendocrine and neuronal derived tissues. In addition, it is internalized after binding to its antigen on SCLC cell lines, making it an excellent reagent for tumor imaging in xenograft models. ab9272 reacts with all three isoforms of NCAM.
经测试应用	适用于: IHC-FrFI, IHC-Fr, Flow Cyt, ICC/IF
种属反应性	与反应: Rat, Human, Zebrafish
免疫原	Tissue, cells or virus corresponding to Human NCAM1. This antibody was raised against a small cell lung carcinoma specimen.
表位	This antibody recognizes an epitope in the NCAM exons 11-13 which is dependent on an intact conformation of the first fibronectin type-III homologous domain encoded by these exons.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	<p>pH: 7.3</p> <p>Preservative: 0.02% Sodium azide</p> <p>Constituents: PBS, 0.1% BSA</p>
纯度	Protein G purified

克隆	单克隆
克隆编号	123C3
骨髓瘤	unknown
同种型	IgG1

### 应用

The Abpromise guarantee      **Abpromise™**承诺保证使用ab9272于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-FrFI		Use at an assay dependent concentration.
IHC-Fr	★☆☆☆☆ (1)	1/500 - 1/1000.
Flow Cyt		Use 1µg for 10 <sup>6</sup> cells. <b>ab170190</b> - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody.
ICC/IF		1/1000.

### 靶标


功能	This protein is a cell adhesion molecule involved in neuron-neuron adhesion, neurite fasciculation, outgrowth of neurites, etc.
序列相似性	Contains 2 fibronectin type-III domains. Contains 5 Ig-like C2-type (immunoglobulin-like) domains.
细胞定位	Secreted and Cell membrane.

### 图片



ab9272 at a dilution of 1/1000, staining NCAM (Alexa 488 secondary at 1/2000) on rat brain tissue (30µm thick coronal sections) in free floating IHC (see protocol link for detailed description). Images showing hippocampal neuron body and processes: [A] 20x objective and [B] punctate cytoplasmic staining; 40x objective. No labeling observed following omission of primary antibody.


Sections were viewed using an Axioplan 2 Imaging microscope (Imaging Associates) fitted with 10x, 20x and 40x Plan-Neofluorobjectives (Zeiss, Germany) and images were taken using a AxioCam Hrm digital camera (Zeiss, Germany) and AxioVision software (Imaging Associates).



Immunohistochemistry (Frozen sections) - Anti-NCAM1 antibody [123C3] (ab9272)

Immunofluorescence staining images of 9 day old zebrafish embryos.


ab9272 reacts with developing neuronal tissue (strong) and muscle. Frozen sample treated with Acetone:Methanol 1:1, antibody diluted 1/100 and incubated for 45 minutes at room temperature.



Immunohistochemistry (Frozen sections) - Anti-NCAM1 antibody [123C3] (ab9272)

Immunofluorescence staining images of 9 day old zebrafish embryos.


ab9272 reacts with developing neuronal tissue (strong) and muscle. Frozen sample treated with Acetone:Methanol 1:1, antibody diluted 1/100 and incubated for 45 minutes at room temperature.



Immunohistochemistry (Frozen sections) - Anti-NCAM1 antibody [123C3] (ab9272)


Immunofluorescence staining images of 9 day old zebrafish embryos.

ab9272 reacts with developing neuronal tissue (strong) and muscle. Frozen sample treated with Acetone:Methanol 1:1, antibody diluted 1/100 and incubated for 45 minutes at room temperature.



Immunocytochemistry/ Immunofluorescence - Anti-NCAM1 antibody [123C3] (ab9272)

Immunostaining of the neural cell adhesion molecule (NCAM) in small cell human lung cancer cell line NCI-H82.



Flow Cytometry - Anti-NCAM1 antibody [123C3] (ab9272)

Overlay histogram showing SH-SY5Y cells stained with ab9272 (red line). The cells were fixed with 4% paraformaldehyde (10 min) and then permeabilized with 0.1% PBS-Tween for 20 min. The cells were then incubated in 1x PBS / 10% normal goat serum / 0.3M glycine to block non-specific protein-protein interactions followed by the antibody (ab9272, 1µg/1x10<sup>6</sup> cells) for 30 min at 22°C. The secondary antibody used was DyLight® 488 goat anti-mouse IgG (H+L) ([ab96879](#)) at 1/500 dilution for 30 min at 22°C. Isotype control antibody (black line) was mouse IgG1 [ICIGG1] ([ab91353](#), 2µg/1x10<sup>6</sup> cells) used under the same conditions. Acquisition of >5,000 events was performed. This antibody gave a positive signal in SH-SY5Y cells fixed with 80% methanol (5 min)/permeabilized with 0.1% PBS-Tween for 20 min used under the same conditions.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors