

Anti-Neuropilin 1 antibody ab25998

★★★★★ [1 Abreviews](#) [14 References](#) [3 图像](#)

概述

| | |
|-------|---|
| 产品名称 | Anti-Neuropilin 1抗体 |
| 描述 | 兔多克隆抗体to Neuropilin 1 |
| 宿主 | Rabbit |
| 特异性 | ab25998 is specific for Neuropilin 1. Immunogen sequence has low homology to Neuropilin 2. |
| 经测试应用 | 适用于: ICC/IF, IHC-P, WB |
| 种属反应性 | 与反应: Human, African green monkey |
| 免疫原 | Synthetic peptide, coupled to carrier protein, corresponding to amino acids within the N-terminal CUB domain of human Neuropilin 1. This sequence is highly conserved in rat and mouse Neuropilin 1. |
| 常规说明 | <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p> |

性能

| | |
|------|---|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |
| 存储溶液 | Preservative: 0.05% Sodium azide Constituents: PBS, 50% Glycerol, 0.1% BSA |
| 纯度 | Immunogen affinity purified |
| 克隆 | 多克隆 |
| 同种型 | IgG |

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab25998于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

| 应用 | Ab评论 | 说明 |
|--------|-----------|---|
| ICC/IF | | 1/200. |
| IHC-P | | 1/100. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol. |
| WB | ★★★★★ (1) | 1/1000. Predicted molecular weight: 102 kDa. |

靶标

| | |
|--------------|--|
| 功能 | <p>The membrane-bound isoform 1 is a receptor involved in the development of the cardiovascular system, in angiogenesis, in the formation of certain neuronal circuits and in organogenesis outside the nervous system. It mediates the chemorepulsant activity of semaphorins. It binds to semaphorin 3A, The PLGF-2 isoform of PGF, The VEGF-165 isoform of VEGF and VEGF-B. Coexpression with KDR results in increased VEGF-165 binding to KDR as well as increased chemotaxis. It may regulate VEGF-induced angiogenesis.</p> <p>The soluble isoform 2 binds VEGF-165 and appears to inhibit its binding to cells. It may also induce apoptosis by sequestering VEGF-165. May bind as well various members of the semaphorin family. Its expression has an averse effect on blood vessel number and integrity.</p> |
| 组织特异性 | <p>The expression of isoforms 1 and 2 does not seem to overlap. Isoform 1 is expressed by the blood vessels of different tissues. In the developing embryo it is found predominantly in the nervous system. In adult tissues, it is highly expressed in heart and placenta; moderately in lung, liver, skeletal muscle, kidney and pancreas; and low in adult brain. Isoform 2 is found in liver hepatocytes, kidney distal and proximal tubules.</p> |
| 序列相似性 | <p>Belongs to the neuropilin family.</p> <p>Contains 2 CUB domains.</p> <p>Contains 2 F5/8 type C domains.</p> <p>Contains 1 MAM domain.</p> |
| 细胞定位 | <p>Secreted and Cell membrane.</p> |


图片



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Neuropilin 1 antibody (ab25998)

IHC image of ab25998 staining in human heart formalin fixed paraffin embedded tissue section, performed on a Leica Bond™ system using the standard protocol F. The section was pre-treated using heat mediated antigen retrieval with sodium citrate buffer (pH6, epitope retrieval solution 1) for 20 mins. The section was then incubated with ab25998, 1µg/ml, for 15 mins at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX.

For other IHC staining systems (automated and non-automated) customers should optimize variable parameters such as antigen retrieval conditions, primary antibody concentration and antibody incubation times.




Western blot - Anti-Neuropilin 1 antibody (ab25998)

All lanes : Anti-Neuropilin 1 antibody (ab25998) at 1/500 dilution

Lane 1 : Untransfected COS-7 cells

Lane 2 : COS-7 cells transfected with HA-tagged Neuropilin 1

Predicted band size: 102 kDa



Immunocytochemistry/ Immunofluorescence - Anti-Neuropilin 1 antibody (ab25998)

ICC/IF image of ab25998 stained SHSY5Y cells. The cells were 4% formaldehyde fixed (10 min) and then incubated in 1%BSA / 10% normal goat serum / 0.3M glycine in 0.1% PBS-Tween for 1h to permeabilise the cells and block non-specific protein-protein interactions. The cells were then incubated with the antibody (ab25998, 5µg/ml) overnight at +4°C. The secondary antibody (green) was Alexa Fluor® 488 goat anti-rabbit IgG (H+L) used at a 1/1000 dilution for 1h. Alexa Fluor® 594 WGA was used to label plasma membranes (red) at a 1/200 dilution for 1h. DAPI was used to stain the cell nuclei (blue) at a concentration of 1.43µM.

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors