

Anti-Neurotrophin 3 antibody ab6203

2 References

概述

产品名称	Anti-Neurotrophin 3抗体
描述	兔多克隆抗体to Neurotrophin 3
宿主	Rabbit
特异性	Less than 1% cross reactivity against NGF, BDNF and NT4/5 by dot blot. Cross reactivity: Known to react with rat, chicken and human NT3. Cross reactivity against other species has not been determined yet.
经测试应用	适用于: IHC-Fr, ELISA, Neutralising, Dot blot, WB
种属反应性	与反应: Human
免疫原	Recombinant full length protein corresponding to Human Neurotrophin 3.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Add glycerol to a final volume of 50% for extra stability and aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 7
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™承诺保证使用ab6203于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-Fr		
ELISA		
Neutralising		
Dot blot		
WB		

应用说明

Immunohistochemistry: 1 ug/ml
Western Blot: 1 ug/ml
Immunoblot: 1 ug/ml
1 site ELISA: 1 ug/ml
Biological Inhibition: 2-10 ug/ml (ED50)

靶标

功能 Seems to promote the survival of visceral and proprioceptive sensory neurons.

组织特异性 Brain and peripheral tissues.

序列相似性 Belongs to the NGF-beta family.

细胞定位 Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors