

Anti-NLRP7 antibody ab117732

★★★★★ [1 Abreviews](#) [1 图像](#)

概述

产品名称	Anti-NLRP7抗体
描述	兔多克隆抗体to NLRP7
宿主	Rabbit
经测试应用	适用于: IHC-P, ELISA, WB
种属反应性	与反应: Human
免疫原	Synthetic peptide conjugated to KLH, corresponding to a region from the N-terminus of Human NLRP7.
阳性对照	Human Spleen tissue
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	<p>pH: 7.4</p> <p>Preservative: 0.02% Sodium azide</p> <p>Constituent: 99% PBS</p>
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab117732于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P		Use a concentration of 5 µg/ml.
ELISA		Use at an assay dependent concentration.
WB	★★★★★ (1)	Use at an assay dependent concentration. Predicted molecular weight: 111 kDa.

靶标

功能	Inhibits CASP1/caspase-1-dependent IL1B secretion.
组织特异性	Expressed in numerous tissues including uterus and ovary, with low levels in heart and brain. Not detected in skeletal muscle.
疾病相关	Defects in NLRP7 are the cause of hydatidiform mole (HYDM) [MIM:231090]. HYDM is an abnormal human pregnancy with no embryo and cystic degeneration of placental villi.
序列相似性	Belongs to the NLRP family. Contains 1 DAPIN domain. Contains 9 LRR (leucine-rich) repeats. Contains 1 NACHT domain.

图片



ab117732 at 5 µg/ml staining NLRP7 in Formalin-Fixed, Paraffin-Embedded Human Spleen.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-NLRP7 antibody (ab117732)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors