

Anti-NUP93 antibody [3332C2a] ab53750

★★★★★ [1 Abreviews](#) [4 References](#) [1 图像](#)

概述

产品名称	Anti-NUP93抗体[3332C2a]
描述	小鼠单克隆抗体[3332C2a] to NUP93
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Recombinant fragment
免疫原	Recombinant Human NUP93 fragment.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at 4°C (stable for up to 12 months). Store at -20°C or -80°C.
存储溶液	pH: 7.40 Preservative: 0.05% Sodium azide Constituents: 1% BSA, PBS
纯度	Protein G purified
克隆	单克隆
克隆编号	3332C2a
同种型	IgG1

应用

**The Abpromise guarantee**      **Abpromise™**承诺保证使用ab53750于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB	★★★★★ (1)	Use at an assay dependent dilution. Predicted molecular weight: 93 kDa.

## 靶标

相关性	Part of the nucleoporin complex; required for correct nuclear pore assembly.
细胞定位	Nucleus, nuclear pore complex. Note: Localizes at the nuclear basket and at or near the nuclear entry to the gated channel of the pore.

## 图片



Anti-NUP93 antibody [3332C2a] (ab53750) + Recombinant immunising protein fragment

**Predicted band size:** 93 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors