

Anti-PAQR8 antibody ab46535

[4 References](#) [2 图像](#)

概述

产品名称	Anti-PAQR8抗体
描述	兔多克隆抗体to PAQR8
宿主	Rabbit
经测试应用	适用于: ICC/IF
种属反应性	与反应: Mouse
免疫原	Synthetic peptide corresponding to PAQR8 (extracellular) conjugated to keyhole limpet haemocyanin.
阳性对照	ICC/IF: p19 cells, NS-1 cells
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.05% Sodium azide Constituent: Whole serum
纯度	Whole antiserum
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab46535于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		Use at an assay dependent concentration.

靶标


相关性	MPRb (Membrane progesterin receptor beta) is a steroid membrane receptor and binds progesterone. MPRb may be involved in oocyte maturation. MPRb is an integral membrane protein and is expressed in brain and testis. It belongs to the ADIPOR family.
细胞定位	Cell Membrane. Multi-pass membrane protein. Note: Colocalizes with a lysosomal protein cathepsin D.

图片



Immunofluorescent analysis of Progesterin Receptor Beta using anti- Progesterin Receptor Beta polyclonal antibody (ab46535) shows staining in NS-1 Cells.

Immunocytochemistry/ Immunofluorescence - Anti-PAQR8 antibody (ab46535)



Immunofluorescent analysis of PAQR8 using anti-PAQR8 polyclonal antibody (ab46535) shows staining in p19 Cells.

Immunocytochemistry/ Immunofluorescence - Anti-PAQR8 antibody (ab46535)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors