

Anti-Paramyxovirus Type 3 antibody ab33637

概述

产品名称	Anti-Paramyxovirus Type 3抗体
描述	鸡多克隆抗体to Paramyxovirus Type 3
宿主	Chicken
特异性	This antibody is monospecific.
经测试应用	适用于: ELISA
种属反应性	与反应: Paramyxovirus sp.
免疫原	Tissue, cells or virus corresponding to Paramyxovirus Type 3.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Add glycerol to a final volume of 50% for extra stability and aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Constituent: Whole serum
纯度	Whole antiserum
纯化说明	This antibody is a sterile filtered serum produced in SPF (Specific Pathogen Free) chicken embryos.
克隆	多克隆
同种型	IgY

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab33637于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		Use at an assay dependent concentration.

靶标

相关性	The family of Paramyxoviridae contains viruses that induce a wide range of distinct clinical illnesses in humans. These include measles virus, which in rare instances is followed by subacute sclerosing panencephalitis (SSPE); mumps virus, which has symptoms of parotitis, orchitis and encephalitis, and the parainfluenza viruses which are respiratory pathogens.
-----	---

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors