

Anti-Parotid secretory protein antibody ab121028

1 图像

概述

产品名称	Anti-Parotid secretory蛋白抗体
描述	山羊多克隆抗体to Parotid secretory蛋白
宿主	Goat
经测试应用	适用于: ICC/IF
种属反应性	与反应: Mouse
免疫原	Synthetic peptide corresponding to Mouse Parotid secretory protein aa 100-200 (internal sequence). Run BLAST with ExpASY Run BLAST with NCBI
阳性对照	ICC/IF: Mouse Submandibular gland cells
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 7.30 Preservative: 0.02% Sodium azide Constituents: 99.48% Tris buffered saline, 0.5% BSA
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用 ab121028 于以下的经测试应用


“应用说明”部分下显示的仅为推荐的起始稀释度；实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		Use a concentration of 4 µg/ml.

靶标

组织特异性	Detected in submandibular gland.
序列相似性	Belongs to the BPI/LBP/Plunc superfamily. Plunc family.
细胞定位	Secreted.

图片



Immunocytochemistry/immunofluorescence analysis of mouse submandibular gland cells labelling parotid secretory protein with ab121028 at 4 µg/mL (blue). Cells were fixed with 4% paraformaldehyde and permeabilized with 0.1% Triton X-100. Cells were counterstained with Ecad at 1/200 dilution (red). Nuclear DNA was labelled with DAPI (green).

Immunocytochemistry/ Immunofluorescence - Anti-Parotid secretory protein antibody (ab121028)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors

