




Anti-PDIA2/PDIP antibody ab110195

1 References **3 图像**

概述

产品名称	Anti-PDIA2/PDIP抗体
描述	山羊多克隆抗体to PDIA2/PDIP
宿主	Goat
经测试应用	适用于: WB, IHC-P
种属反应性	与反应: Human 预测可用于: Mouse, Rat 
免疫原	Synthetic peptide corresponding to Human PDIA2/PDIP aa 300-400 (internal sequence). Database link: Q13087 <div>  Run BLAST with  Run BLAST with </div>
阳性对照	Human Heart and Human Small Intestine tissue. Human Pancreas lysate
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	pH: 7.30 Preservative: 0.02% Sodium azide Constituents: 0.5% BSA, Tris buffered saline
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用 ab110195 于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度; 实际最佳的稀释度/浓度应由使用者检定。

应用	Ab 评论	说明
WB		Use a concentration of 0.01 - 0.03 µg/ml. Predicted molecular weight: 57 kDa.
IHC-P		Use a concentration of 3.75 µg/ml.

靶标


功能	Acts as an intracellular estrogen-binding protein. May be involved in modulating cellular levels and biological functions of estrogens in the pancreas. May act as a chaperone that inhibits aggregation of misfolded proteins.
组织特异性	Highly expressed in pancreas (at protein level).
序列相似性	Belongs to the protein disulfide isomerase family. Contains 2 thioredoxin domains.
翻译后修饰	The disulfide-linked homodimer exhibits an enhanced chaperone activity. Glycosylated.
细胞定位	Endoplasmic reticulum lumen.

图片




ab110195 at 3.75 µg/ml staining PDIA2/PDIP in Human Heart
Formalin-Fixed, Paraffin-Embedded tissue.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-PDIA2/PDIP antibody (ab110195)



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-PDIA2/PDIP antibody (ab110195)

ab110195 at 3.75 µg/ml staining PDIA2/PDIP in Human small intestine Formalin-Fixed, Paraffin-Embedded tissue.



Western blot - Anti-PDIA2/PDIP antibody (ab110195)

Anti-PDIA2/PDIP antibody (ab110195) at 0.03 µg/ml + Human Pancreas lysate at 35 µg

Developed using the ECL technique.

Predicted band size: 57 kDa

Primary incubation was 1 hour

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors