

PE Anti-CD1d antibody [51.1] ab95733

[1 References](#) [1 图像](#)

概述

产品名称	PE Anti-CD1d抗体[51.1]
描述	PE小鼠单克隆抗体[51.1] to CD1d
宿主	Mouse
偶联物	PE. Ex: 488nm, Em: 575nm
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Human
免疫原	Fc-fusion protein of Human CD1d expressed in Mouse NSO cells
阳性对照	Normal Human peripheral blood cells
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	<p>pH: 7.20</p> <p>Preservative: 0.09% Sodium azide</p> <p>Constituents: 0.1% Gelatin, PBS</p>
纯度	Containing horse serum
克隆	Protein G purified
克隆编号	单克隆
同种型	51.1
	IgG2b

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab95733于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use 5µl for 10 ⁶ cells. ab91532 - Mouse monoclonal IgG2b, is suitable for use as an isotype control with this antibody.

靶标

功能	Antigen-presenting protein that binds self and non-self glycolipids and presents them to T-cell receptors on natural killer T-cells.
组织特异性	Expressed on cortical thymocytes, on certain T-cell leukemias, and in various other tissues.
序列相似性	Contains 1 Ig-like (immunoglobulin-like) domain.
细胞定位	Cell membrane. Endosome membrane. Lysosome membrane. Subject to intracellular trafficking between the cell membrane, endosomes and lysosomes.

图片



Staining of normal Human peripheral blood cells with Anti-Human CD19 FITC and Mouse IgG2b ? Isotype Control PE (left) or ab95733 (right). Cells in the lymphocyte gate were used for analysis.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors