

PE Anti-CD200 / OX2 antibody [MRC OX90] ab33739

[1 References](#) [1 图像](#)

概述

产品名称	PE Anti-CD200 / OX2抗体[MRC OX90]
描述	PE大鼠单克隆抗体[MRC OX90] to CD200 / OX2
宿主	Rat
偶联物	PE. Ex: 488nm, Em: 575nm
特异性	Recognises the mouse CD200 cell surface antigen, also known as OX2.
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Mouse
免疫原	Fusion protein corresponding to Mouse CD200/ OX2. Database link: O54901
阳性对照	CD200 is expressed by splenic B lymphocytes, follicular dendritic cells, splenic endothelium and by neurons.
常规说明	Purified IgG conjugated to R. Phycoerythrin (RPE). The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing. If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.09% Sodium azide Constituents: PBS, 1% BSA
纯度	Protein G purified
纯化说明	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant.
克隆	单克隆
克隆编号	MRC OX90

骨髓瘤 Y3/Ag1.2.3

同种型 IgG2a

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab33739于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use 10µl for 10 ⁶ cells. Use 10µl for 10 ⁶ cells in 100µl. The Fc region of monoclonal antibodies may bind non-specifically to cells expressing low affinity fc receptors.


靶标

功能 Costimulates T-cell proliferation. May regulate myeloid cell activity in a variety of tissues.

序列相似性 Contains 1 Ig-like C2-type (immunoglobulin-like) domain.
Contains 1 Ig-like V-type (immunoglobulin-like) domain.

细胞定位 Cell membrane.

图片



ab33739 staining of mouse splenocytes.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors