

PE Anti-CD28 antibody [YTH 913.12] ab30387

1 图像

概述

产品名称	PE Anti-CD28抗体[YTH 913.12]
描述	PE大鼠单克隆抗体[YTH 913.12] to CD28
宿主	Rat
偶联物	PE. Ex: 488nm, Em: 575nm
特异性	It has been reported that this clone also recognises an epitope of CD28 expressed by NK cells, which is not recognised by clones of 9.3 or CD28.2.
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Human
免疫原	Tissue, cells or virus corresponding to Human CD28. Human peripheral blood T cells.
阳性对照	Human peripheral blood lymphocytes.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	<p>pH: 7.40</p> <p>Preservative: 0.09% Sodium azide</p> <p>Constituents: PBS, 1% BSA</p>
纯度	Protein G purified
纯化说明	Purified from tissue culture supernatant.
克隆	单克隆
克隆编号	YTH 913.12
骨髓瘤	Y3/Ag1.2.3

应用

The Abpromise guarantee

Abpromise™承诺保证使用ab30387于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		

应用说明

Flow Cyt: use neat.
Use 10ul of the suggested working dilution to label 106 cells or 100ul human whole blood.

Not yet tested in other applications.
Optimal dilutions/concentrations should be determined by the end user.

靶标

功能

Involved in T-cell activation, the induction of cell proliferation and cytokine production and promotion of T-cell survival.

组织特异性

Expressed in T-cells and plasma cells, but not in less mature B-cells.


序列相似性

Contains 1 Ig-like V-type (immunoglobulin-like) domain.

细胞定位

Membrane.

图片



Profile of human peripheral blood lymphocytes stained with rat anti-human CD28 [YTH 913.12] phycoerythrin conjugate antibody.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors