

PE Anti-CD30 antibody [2SH12-5F] ab25187

1 图像

概述

产品名称	PE Anti-CD30抗体[2SH12-5F]
描述	PE亚美尼亚仓鼠单克隆抗体[2SH12-5F] to CD30
宿主	Armenian hamster
偶联物	PE. Ex: 488nm, Em: 575nm
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Mouse
免疫原	Recombinant CD30 extracellular domain-mouse IgG1 fusion protein (Mouse)
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C. Store In the Dark.
存储溶液	<p>pH: 7.3</p> <p>Preservative: 0.09% Sodium azide</p> <p>Constituents: PBS, 16% Sucrose</p> <p>Stabilising agent.</p>
纯度	Affinity purified
克隆	单克隆
克隆编号	2SH12-5F
同种型	IgG
轻链类型	kappa

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab25187于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use 0.2µg for 10 ⁶ cells.

靶标

功能	Receptor for TNFSF8/CD30L. May play a role in the regulation of cellular growth and transformation of activated lymphoblasts. Regulates gene expression through activation of NF-kappa-B.
序列相似性	Contains 6 TNFR-Cys repeats.
翻译后修饰	Phosphorylated on serine and tyrosine residues.
细胞定位	Cytoplasm and Cell membrane.

图片



Flow cytometry analysis staining CD30 in NIH/3T3 and mCD30-Ig transfected 3T3 cells using ab25187 at a dilution of 0.1 µg/10⁶ cells. Hamster IgG-PE was used as an isotype control.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors