

### PE Anti-CD32B + CD32A antibody [AT10] ab30357

[8 References](#) [1 图像](#)

#### 概述

产品名称	PE Anti-CD32B + CD32A抗体[AT10]
描述	PE小鼠单克隆抗体[AT10] to CD32B + CD32A
宿主	Mouse
偶联物	PE. Ex: 488nm, Em: 575nm
特异性	According to the CD32 Cluster workshop (Leucocyte Typing V, pages 823 and 828), clone AT10 reacts equally well with both Fc gamma RII a and b-expressing cells.
经测试应用	<b>适用于:</b> Flow Cyt
种属反应性	<b>与反应:</b> Human
免疫原	Tissue, cells or virus corresponding to Human CD32A. K562 cell line
阳性对照	Human peripheral blood granulocytes
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.20 Preservative: 0.09% Sodium azide Constituents: PBS, 1% BSA
纯度	Protein A purified
克隆	单克隆
克隆编号	AT10
骨髓瘤	NS1
同种型	IgG1



If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors