

PE Anti-CD90 / Thy1 antibody [G7] ab24904

★★★★★ [1 Abreviews](#) [8 References](#) [1 图像](#)

概述

产品名称	PE Anti-CD90 / Thy1抗体[G7]
描述	PE大鼠单克隆抗体[G7] to CD90 / Thy1
宿主	Rat
偶联物	PE. Ex: 488nm, Em: 575nm
特异性	ab24904 is specific for Thy 1.1 and Thy 1.2 alloantigens (Thy 1 epitope region A).
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Mouse
免疫原	Tissue, cells or virus corresponding to Mouse CD90/ Thy1. Mouse T cell hybridoma C6/G8
表位	Thy 1 epitope region A
常规说明	<p>ab24904 stimulates T-cell proliferation and IL2 secretion, via signalling through the T-cell receptor/CD3 complex. ab24904 has also been reported to promote apoptosis of thymocytes and CTL clones and to mediate adhesion of thymocytes to thymic stroma.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C. Store In the Dark.
存储溶液	<p>pH: 7.3</p> <p>Preservative: 0.09% Sodium azide</p> <p>Constituents: PBS, 16% Sucrose</p> <p>Stabilising agent.</p>

纯度	Affinity purified
Primary antibody说明	ab24904 stimulates T-cell proliferation and IL2 secretion, via signalling through the T-cell receptor/CD3 complex. ab24904 has also been reported to promote apoptosis of thymocytes and CTL clones and to mediate adhesion of thymocytes to thymic stroma.
克隆	单克隆
克隆编号	G7
同种型	IgG2c
轻链类型	kappa

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab24904于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt	★★★★★ (1)	Use 2µg for 10 ⁶ cells. ab154450 - Rat monoclonal IgG2c, is suitable for use as an isotype control with this antibody.

靶标

功能	May play a role in cell-cell or cell-ligand interactions during synaptogenesis and other events in the brain.
序列相似性	Contains 1 Ig-like V-type (immunoglobulin-like) domain.
细胞定位	Cell membrane.

图片



Flow cytometry analysis staining CD90 / Thy1 in mouse splenocytes using ab24904 at a dilution of 0.03 µg/10⁶ cells.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors