

PE Anti-HLA Class I antibody [W6/32] ab33257

8 References

概述

产品名称	PE Anti-HLA Class I 抗体[W6/32]
描述	PE小鼠单克隆抗体[W6/32] to HLA Class I
宿主	Mouse
偶联物	PE. Ex: 488nm, Em: 575nm
特异性	This antibody recognises an antigenic determinant shared among products of the HLA A, B and C loci.
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Cow, Cat, Human, Monkey, Baboon, Cynomolgus monkey, Macaque monkey
免疫原	Tissue, cells or virus corresponding to Human HLA Class I. Purified human tonsil lymphocyte membranes.
阳性对照	Human peripheral blood lymphocytes.
常规说明	<p>50 tests will be supplied typically at a volume of 0.5ml at a concentration of 0.1mg/ml (concentration prior to conjugation).</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.09% Sodium azide Constituents: PBS, 1% BSA, 5% Sucrose
纯度	Protein G purified
纯化说明	This antibody was purified from tissue culture supernatant.
克隆	单克隆

克隆编号	W6/32
骨髓瘤	NS1/1-Ag4-1
同种型	IgG2a

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab33257于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		<p>Use 10µl for 10⁶ cells.</p> <p>ab91363 - Mouse monoclonal IgG2a, is suitable for use as an isotype control with this antibody.</p> <p>Use 10µl to label 10⁶ cells or 100µl whole blood.</p>

靶标

相关性	HLA Class I is involved in the presentation of foreign antigens to the immune system.
细胞定位	Plasma membrane

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors