

PE Mouse IgG1 monoclonal - Isotype control ab81200

2 References

概述

产品名称	PE小鼠IgG1单克隆抗体-同型对照
偶联物	PE. Ex: 488nm, Em: 575nm
共轭说明	R-PE excitation wavelength 488 nm, emission wavelength 575 nm is the correct sentence.
经测试应用	适用于: Flow Cyt
常规说明	<p>Reagents should be brought to room temperature ($22 \pm 3^{\circ}\text{C}$) before use.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C . Store at $+4^{\circ}\text{C}$ short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
存储溶液	<p>pH: 7.3</p> <p>Preservative: 0.1% Sodium azide</p> <p>Constituents: 0.4% BSA, 1.74% Sodium chloride, 0.328% Sodium phosphate</p>
纯化说明	Purified from mouse myeloma using $0.22\ \mu\text{m}$ filtration
克隆	单克隆
克隆编号	TIB9
同种型	IgG1

应用

The Abpromise guarantee **Abpromise™** 承诺保证使用ab81200于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use at an assay dependent concentration. Nonspecific binding in flow cytometric methods of various diseases.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors