

PE Anti-NCAM1 antibody [MEM-188] ab18277

3 References

概述

产品名称	PE Anti-NCAM1抗体[MEM-188]
描述	PE小鼠单克隆抗体[MEM-188] to NCAM1
宿主	Mouse
偶联物	PE. Ex: 488nm, Em: 575nm
特异性	Human NCAM. Reacts with the 180 kDa isoform characteristic for leukocytes. Reactivity with other NCAM isoforms has not been tested.
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Mouse, Human, Monkey
免疫原	Tissue, cells or virus corresponding to NCAM1. KG-1 human acute myelogenous leukaemia cell line
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C. Store In the Dark.
存储溶液	<p>pH: 7.40</p> <p>Preservative: 0.097% Sodium azide</p> <p>Constituents: PBS, BSA</p>
纯化说明	Purified by size exclusion chromatography.
克隆	单克隆
克隆编号	MEM-188
同种型	IgG2a

应用

The Abpromise guarantee Abpromise™承诺保证使用ab18277于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use at an assay dependent concentration. Designed for analysis of human blood cells using 20ul reagent per 100ul whole blood or 10 ⁶ cells in a suspension. ab91363 - Mouse monoclonal IgG2a, is suitable for use as an isotype control with this antibody.

靶标

功能	This protein is a cell adhesion molecule involved in neuron-neuron adhesion, neurite fasciculation, outgrowth of neurites, etc.
序列相似性	Contains 2 fibronectin type-III domains. Contains 5 Ig-like C2-type (immunoglobulin-like) domains.
细胞定位	Secreted and Cell membrane.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors