

### PE/Cy5® Anti-CD90 / Thy1 antibody [5E10] ab95698

[7 References](#) [1 图像](#)

#### 概述

<b>产品名称</b>	PE/Cy5® Anti-CD90 / Thy1 抗体[5E10]
<b>描述</b>	PE/Cy5®小鼠单克隆抗体[5E10] to CD90 / Thy1
<b>宿主</b>	Mouse
<b>偶联物</b>	PE/Cy5®. Ex: 496nm, Em: 670nm
<b>经测试应用</b>	<b>适用于:</b> Flow Cyt
<b>种属反应性</b>	<b>与反应:</b> Human
<b>免疫原</b>	Tissue, cells or virus corresponding to Human CD90/ Thy1. HEL cells
<b>阳性对照</b>	Flow cyt: HEL cells.
<b>常规说明</b>	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 性能

<b>形式</b>	Liquid
<b>存放说明</b>	Shipped at 4°C. Store at +4°C.
<b>存储溶液</b>	pH: 7.20 Preservative: 0.09% Sodium azide Constituents: BSA, PBS
<b>纯度</b>	Protein G purified
<b>克隆</b>	单克隆
<b>克隆编号</b>	5E10
<b>同种型</b>	IgG1
<b>轻链类型</b>	kappa

## 应用

**The Abpromise guarantee** **Abpromise™**承诺保证使用ab95698于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use at an assay dependent concentration. 5 µL (0.25 µg)/test


## 靶标

**功能** May play a role in cell-cell or cell-ligand interactions during synaptogenesis and other events in the brain.

**序列相似性** Contains 1 Ig-like V-type (immunoglobulin-like) domain.

**细胞定位** Cell membrane.

## 图片



Flow cytometry analysis of HEL (human bone marrow erythroleukemia cell line) cells labelling CD90/Thy1 with Mouse IgG1 K Isotype Control PE-Cyanine5 (blue histogram) or **ab35698** (purple histogram). Total viable cells were used for analysis.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors