

PE/Cy7® Anti-CD3 antibody [UCHT1] ab81992

1 图像

概述

产品名称	PE/Cy7® Anti-CD3抗体[UCHT1]
描述	PE/Cy7®小鼠单克隆抗体[UCHT1] to CD3
宿主	Mouse
偶联物	PE/Cy7®. Ex: 496nm, Em: 774nm
特异性	Reacts with the 20 kDa e chain of the CD3/TCR complex
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Human
免疫原	Tissue, cells or virus corresponding to Human CD3. Human infant thymocytes and peripheral blood lymphocytes from a Sézary Syndrome donor
阳性对照	Flow Cyt: Human peripheral blood lymphocytes.
常规说明	<p>This product or portions thereof is manufactured under license from Carnegie Mellon University under U.S. Patent Number 5, 268, 486 and related patents. Cy® and CyDye® are trademarks of Cytiva.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C. Do Not Freeze. Store In the Dark.
存储溶液	<p>pH: 7.20</p> <p>Preservative: 0.09% Sodium azide</p> <p>Constituents: 0.1% Gelatin, 0.88% Sodium chloride, 0.12% Monobasic dihydrogen sodium phosphate</p>
纯度	Affinity purified

纯化说明	Purified from TCS
克隆	单克隆
克隆编号	UCHT1
同种型	IgG1


应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab81992于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use 5µl for 10 ⁵⁻⁸ cells. Use 5µl in a final volume of 100µl. ab154457 - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody.

图片



Human peripheral blood lymphocytes stained with 5µL (0.5µg) ab81992 (solid line) or 0.5µg PE-Cy7-conjugated mouse IgG1 isotype control (dashed line).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors