


### PerCP Anti-CD28 antibody [CD28.2] ab88034

#### 7 References

#### 概述

产品名称	PerCP荧光Anti-CD28抗体[CD28.2]
描述	PerCP荧光小鼠单克隆抗体[CD28.2] to CD28
宿主	Mouse
偶联物	PerCP. Ex: 482nm, Em: 675nm
特异性	Reacts with CD28, a disulfide-linked homodimeric type I glycoprotein (monomer of Mw 44 kDa) which is a critical costimulatory receptor of T cells.
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Human 预测可用于: Non human primates 
免疫原	Tissue, cells or virus corresponding to CD28. DC28.1.3.3 hybridoma cells
阳性对照	Human blood cells
常规说明	<p>The purified antibody (&gt;95% by SDS-PAGE) is conjugated with Peridinin-chlorophyll-protein complex (PerCP) under optimum conditions.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	<p>pH: 7.4</p> <p>Preservative: 0.097% Sodium azide</p> <p>Constituent: 0.2% BSA</p> <p>0.2% BSA (high-grade, protease free)</p>

纯度	Size exclusion
纯化说明	The conjugate is purified by size-exclusion chromatography and adjusted for direct use.
克隆	单克隆
克隆编号	CD28.2
同种型	IgG1

## 应用

**The Abpromise guarantee** Abpromise™承诺保证使用ab88034于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use 10µl for 10 <sup>6</sup> cells.

## 靶标

功能	Involved in T-cell activation, the induction of cell proliferation and cytokine production and promotion of T-cell survival.
组织特异性	Expressed in T-cells and plasma cells, but not in less mature B-cells.
序列相似性	Contains 1 Ig-like V-type (immunoglobulin-like) domain.
细胞定位	Membrane.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors