


Anti-Perilipin 2 antibody ab52356

★★★★★ [3 Abreviews](#) [31 References](#) [3 图像](#)

概述

| | |
|-------|---|
| 产品名称 | Anti-Perilipin 2抗体 |
| 描述 | 兔多克隆抗体to Perilipin-2 |
| 宿主 | Rabbit |
| 经测试应用 | 适用于: ICC, WB, IHC-P |
| 种属反应性 | 与反应: Mouse, Human 预测可用于: Rat  |
| 免疫原 | Synthetic peptide corresponding to Mouse Perilipin-2 (C terminal). |
| 常规说明 | <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p> |

性能

| | |
|------|---|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |
| 存储溶液 | pH: 7.40 Constituents: 30% Glycerol, PBS |
| 纯度 | Immunogen affinity purified |
| 克隆 | 多克隆 |
| 同种型 | IgG |

应用

The Abpromise guarantee **Abpromise™** 承诺保证使用ab52356于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

| 应用 | Ab评论 | 说明 |
|-------|-----------|--|
| ICC | | Use at an assay dependent concentration. |
| WB | ★★★★★ (2) | Use a concentration of 2 µg/ml. Detects a band of approximately 48 kDa (predicted molecular weight: 48 kDa). |
| IHC-P | | 1/200. |

靶标


| | |
|-------|---|
| 功能 | May be involved in development and maintenance of adipose tissue. |
| 组织特异性 | Milk lipid globules. |
| 序列相似性 | Belongs to the perilipin family. |
| 翻译后修饰 | Acylated; primarily with C14, C16 and C18 fatty acids. |
| 细胞定位 | Membrane. |

图片




Anti-Perilipin 2 antibody (ab52356) at 2 µg/ml + mouse liver lysate

Predicted band size: 48 kDa
Observed band size: 48 kDa



Immunohistochemical analysis of mouse liver tissue staining Perilipin-2 using ab52356.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Perilipin 2 antibody (ab52356)



Ab52356 staining Perilipin-2 in HeLa cells at 1/50 dilution. Hoechst 33258 was used as a nuclear counterstain.

Immunocytochemistry - Anti-Perilipin 2 antibody (ab52356)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors