


Anti-Periostin antibody ab92460

★★★★★ 1 Abreviews 18 References 1 图像

|       |   |
|-------|---|
| 概述    |   |
| 产品名称  | Anti-Periostin抗体  |
| 描述    | 兔多克隆抗体to Periostin  |
| 宿主    | Rabbit  |
| 经测试应用 | 适用于: WB, ICC/IF, IHC-P  |
| 种属反应性 | 与反应: Mouse, Rat, Chicken, Human<br>预测可用于: Bird, Fish, Mammals, Amphibian   |
| 免疫原   | Synthetic peptide corresponding to Mouse Periostin. Synthetic peptide from the fasciclin domain 1 of mouse periostin.<br>Database link: <a href="#">Q62009</a>  |
| 常规说明  | <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p> |
| 性能    |   |
| 形式    | Liquid  |
| 存放说明  | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.   |
| 存储溶液  | Constituent: PBS  |
| 纯度    | Immunogen affinity purified   |
| 克隆    | 多克隆   |
| 同种型   | IgG   |
| 应用    |   |

The Abpromise guarantee      **Abpromise™**承诺保证使用ab92460于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

| 应用     | Ab评论      | 说明  |
|--------|-----------|---|
| WB     | ★★★★★ (1) | 1/1000. Predicted molecular weight: 93 kDa. |
| ICC/IF |           | Use at an assay dependent concentration.    |
| IHC-P  |           | Use at an assay dependent concentration.    |

## 靶标

|       |   |
|-------|---|
| 功能    | Binds to heparin. Induces cell attachment and spreading and plays a role in cell adhesion. May play a role in extracellular matrix mineralization.  |
| 组织特异性 | Widely expressed with highest levels in aorta, stomach, lower gastrointestinal tract, placenta, uterus and breast. Up-regulated in epithelial ovarian tumors. Not expressed in normal ovaries. Also highly expressed at the tumor periphery of lung carcinoma tissue but not within the tumor. Overexpressed in breast cancers. |
| 序列相似性 | Contains 1 EMI domain.<br>Contains 4 FAS1 domains.  |
| 翻译后修饰 | Gamma-carboxyglutamate residues are formed by vitamin K dependent carboxylation. These residues are essential for the binding of calcium.   |
| 细胞定位  | Secreted > extracellular space > extracellular matrix.  |

## 图片



Anti-Periostin antibody (ab92460) at 1/1000 dilution + rat lung lysate

**Predicted band size:** 93 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors