

Anti-Peripherin antibody ab106276

[2 References](#) [1 图像](#)

概述

产品名称	Anti-Peripherin抗体
描述	鸡多克隆抗体to Peripherin
宿主	Chicken
经测试应用	适用于: ICC/IF
种属反应性	与反应: Rat
免疫原	Rat peripherin was expressed in bacteria and purified from inclusion bodies.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.065% Sodium azide Constituents: PBS, Whole serum
纯度	Whole antiserum
克隆	多克隆
同种型	IgY

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab106276于以下的经测试应用


“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		1/1000.

靶标

功能	Class-III neuronal intermediate filament protein.
序列相似性	Belongs to the intermediate filament family.

图片



Immunocytochemistry/ Immunofluorescence - Anti-Peripherin antibody (ab106276)

ab106276, at 1/2000, staining Peripherin (green channel) and Neurofilament alpha-interneixin/NF66 (green channel) in Rat mixed neuron/glia cultures by Immunofluorescence. These cultures contain mostly neurons which are rich in alpha interneixin, and a subgroup which have a large amount of peripherin also, such as the prominent cell in the middle of the micrograph. Since this cell expresses large amounts of peripherin and alpha interneixin, the green and red signals superimpose to produce a golden cell. Blue is a DNA stain.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors