

# Anti-Placental alkaline phosphatase (PLAP) antibody [H7E8] ab118868

### 1 图像

#### 概述

产品名称	Anti-Placental alkaline phosphatase (PLAP)抗体[H7E8]
描述	小鼠单克隆抗体[H7E8] to Placental alkaline phosphatase (PLAP)
宿主	Mouse
经测试应用	适用于: IHC-P
种属反应性	与反应: Human
免疫原	Purified full length native Human Placental alkaline phosphatase (PLAP)
阳性对照	Human brain, cerebellum, tissue.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at -20°C or -80°C. Avoid freeze / thaw cycle. Store undiluted.
存储溶液	pH: 7.40 Preservative: 0.02% Sodium azide Constituent: 99% PBS
纯度	Protein A purified
克隆	单克隆
克隆编号	H7E8
同种型	IgG2a

#### 应用

## The Abpromise guarantee

**Abpromise™** 承诺保证使用 ab118868 于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度；实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P		Use a concentration of 10 µg/ml. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

## 靶标


### 序列相似性

Belongs to the alkaline phosphatase family.

### 细胞定位

Cell membrane.

## 图片



ab118868 at 10 µg/ml staining Placental alkaline phosphatase (PLAP) in Paraffin-embedded Human cerebellum tissue by immunohistochemistry. After incubation with the primary antibody, slides were incubated with biotinylated secondary antibody, followed by alkaline phosphatase-streptavidin and chromogen.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Placental alkaline phosphatase (PLAP) antibody [H7E8] (ab118868)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors