

Anti-Pma1 antibody [40B7] ab4645

★★★★★ **1 Abreviews** **28 References**

概述

产品名称	Anti-Pma1抗体[40B7]
描述	小鼠单克隆抗体[40B7] to Pma1
宿主	Mouse
特异性	Reacts with Pma1p based on the size of the band detected on western blots and the immunofluorescence localization pattern.
经测试应用	适用于: ICC/IF, WB
种属反应性	与反应: <i>Saccharomyces cerevisiae</i>
免疫原	Tissue, cells or virus corresponding to Pma1. Yeast nuclear prep
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.065% Sodium azide Constituent: Tissue culture supernatant
纯度	Tissue culture supernatant
纯化说明	Sterile filtered.
克隆	单克隆
克隆编号	40B7
同种型	IgG

The Abpromise guarantee **Abpromise™**承诺保证使用ab4645于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		Use at an assay dependent concentration.
WB	★★★★★ (1)	Use at an assay dependent concentration. (cell lysates, ECL detection). For non-ECL western detection methods, 1/1000 - 1/5000. Also detects Pma1p-containing fusion proteins on western blots.

靶标

相关性 Pma1p is a plasma membrane ATPase.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors