

Anti-POLR2F antibody [2780C3a] ab51362

1 图像

概述	
产品名称	Anti-POLR2F 抗体[2780C3a]
描述	小鼠单克隆抗体[2780C3a] to POLR2F
宿主	Mouse
经测试应用	适用于: WB
免疫原	Recombinant fragment (Human)
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>
性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
存储溶液	pH: 7.40 Preservative: 0.05% Sodium azide Constituents: 1% BSA, PBS
纯度	Protein G purified
纯化说明	Filtered through a 0.22 µm membrane.
克隆	单克隆
克隆编号	2780C3a
同种型	IgG1
应用	


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		Use at an assay dependent concentration. Predicted molecular weight: 14 kDa.

靶标

相关性	POLR2F belongs to the archaeal rpoK/eukaryotic RPB6 RNA polymerase subunit family and is the sixth largest subunit of RNA polymerase II. This polymerase, along with two other DNA-directed RNA polymerases, is responsible for synthesizing messenger RNA in eukaryotes. Three distinct zinc-containing RNA polymerases are found in eukaryotic nuclei: polymerase I for the ribosomal RNA precursor, polymerase II for the mRNA precursor, and polymerase III for 5S and tRNA genes. Each class of RNA polymerase is assembled from 9 to 15 different polypeptides. In yeast, this polymerase subunit, in combination with at least two other subunits, forms a structure that stabilizes the transcribing polymerase on the DNA template.
细胞定位	Nuclear

图片



Anti-POLR2F antibody [2780C3a] (ab51362) + immunizing peptide

Predicted band size: 14 kDa

Observed band size: 38 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors