

### Anti-PON3 antibody [EPR2902(2)] ab108964

重组 RabMAb

#### 2 图像

#### 概述

产品名称	Anti-PON3抗体[EPR2902(2)]
描述	兔单克隆抗体[EPR2902(2)] to PON3
宿主	Rabbit
经测试应用	<b>适用于:</b> WB <b>不适用于:</b> Flow Cyt, ICC/IF, IHC-P or IP
种属反应性	<b>与反应:</b> Mouse, Human
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
阳性对照	Human fetal liver, A549, HepG2, Mouse kidney, and Mouse liver lysates
常规说明	<p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none"><li>- High batch-to-batch consistency and reproducibility</li><li>- Improved sensitivity and specificity</li><li>- Long-term security of supply</li><li>- Animal-free production</li></ul> <p>For more information <a href="#">see here</a>.</p> <p>Our RabMAb<sup>®</sup> technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to <a href="#">RabMAb<sup>®</sup> patents</a>.</p> <p>Rat: We have preliminary internal testing data to indicate this antibody may not react with this species. Please contact us for more information.</p>

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
存储溶液	pH: 7.20 Preservative: 0.05% Sodium azide Constituents: 0.1% BSA, 40% Glycerol (glycerin, glycerine), 9.85% Tris glycine, 50% Tissue culture supernatant
纯度	Tissue culture supernatant
克隆	单克隆

克隆编号 EPR2902(2)

同种型 IgG

## 应用

**The Abpromise guarantee** **Abpromise™**承诺保证使用ab108964于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/1000 - 1/10000. Detects a band of approximately 40 kDa (predicted molecular weight: 40 kDa).

应用说明 Is unsuitable for Flow Cyt, ICC/IF, IHC-P or IP.

## 靶标


**功能** Has low activity towards the organophosphate paraoxon and aromatic carboxylic acid esters. Rapidly hydrolyzes lactones such as statin prodrugs (e.g. lovastatin). Hydrolyzes aromatic lactones and 5- or 6-member ring lactones with aliphatic substituents but not simple lactones or those with polar substituents.

**序列相似性** Belongs to the paraoxonase family.

**翻译后修饰** The signal sequence is not cleaved.

**细胞定位** Secreted > extracellular space.

## 图片



**All lanes** : Anti-PON3 antibody [EPR2902(2)] (ab108964) at 1/1000 dilution

**Lane 1** : Human fetal liver lysate

**Lane 2** : A549 lysate

**Lane 3** : HepG2 lysate


**Lane 4** : Mouse kidney lysate

**Lane 5** : Mouse liver lysate

Lysates/proteins at 10 µg per lane.

**Predicted band size:** 40 kDa

Why choose a recombinant antibody?



- Research with confidence**  
Consistent and reproducible results
- Long-term and scalable supply**  
Recombinant technology
- Success from the first experiment**  
Confirmed specificity
- Ethical standards compliant**  
Animal-free production

Anti-PON3 antibody [EPR2902(2)] (ab108964)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors