

Anti-Porcine endothelial cells antibody [MIL11] ab34023

1 References

概述

产品名称	Anti-Porcine endothelial cells抗体[MIL11]
描述	小鼠单克隆抗体[MIL11] to Porcine endothelial cells
宿主	Mouse
特异性	ab34023 recognises an antigen expressed by most endothelial cells in the pig. Expression is strongest upon intestinal endothelial cells, and staining is also seen upon endothelial cells in the lung, colon, trachea, thymus, lymph nodes, kidney and skin.
经测试应用	适用于: IHC-Fr
种属反应性	与反应: Pig
免疫原	Tissue, cells or virus corresponding to Pig Porcine endothelial cells. Plastic adherent porcine peripheral blood lymphocytes.
阳性对照	Porcine endothelial cell lysate
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
存储溶液	Preservative: 0.09% Sodium azide Constituent: PBS
纯度	Tissue culture supernatant
克隆	单克隆
克隆编号	MIL11
同种型	IgE

应用

The Abpromise guarantee Abpromise™承诺保证使用ab34023于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-Fr		1/20 - 1/100.

靶标

相关性 This antibody recognises an antigen expressed by most endothelial cells in the pig. Expression is strongest upon intestinal endothelial cells, and staining is also seen upon endothelial cells in the lung, colon, trachea, thymus, lymph nodes, kidney and skin. In blood vessels, venous endothelial cells but not arterial endothelial cell staining is observed.

细胞定位 Cell Membrane

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors