

Anti-PPAN antibody ab122844

1 References **3 图像**


概述

产品名称	Anti-PPAN抗体
描述	兔多克隆抗体to PPAN
宿主	Rabbit
经测试应用	适用于: IHC-P, WB, ICC/IF
种属反应性	与反应: Human
免疫原	antigen sequence: RWEMDRGRGR LCDQKFPKTK DKSQGAQARR GPRGASRDGG, corresponding to amino acids 421-460 of Human PPAN. Run BLAST with Expasy Run BLAST with NCBI
阳性对照	RT4, U-251 MG and Human Tonsil lysates, Human oral mucosa tissue, U-251 MG cells
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能


形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol (glycerin, glycerine)
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用



ab122844 at 1/200 dilution staining PPAN in Paraffin-embedded Human oral mucosa tissue by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-PPAN antibody (ab122844)



Immunofluorescent staining of Human cell line U-251 MG using ab122844 at a concentration of 4 µg/ml. ab122844 staining is shown in green.

Immunocytochemistry/ Immunofluorescence - Anti-PPAN antibody (ab122844)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors