

Anti-Proteasome 19S S4/PSMC1 antibody ab3317

★★★★★ [1 Abreviews](#) [4 References](#) [1 图像](#)

概述

产品名称	Anti-Proteasome 19S S4/PSMC1抗体
描述	兔多克隆抗体to Proteasome 19S S4/PSMC1
宿主	Rabbit
特异性	Detects proteasome 19S subunit S4.
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Full length protein corresponding to Proteasome 19S S4/PSMC1. Yeast proteasome 19S subunit S4.
阳性对照	WB: HeLa cell lysate.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Constituents: 0.1% BSA, 99% PBS
纯度	IgG fraction
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab3317于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/1000. Detects a band of approximately 56 kDa (predicted molecular weight: 49 kDa). This antibody detects a 56 kDa protein representing proteasome 19S subunit S4 in HeLa cell lysate. This antibody detects, to a lesser extent, an ~35 kDa protein which could correspond to subunit degradation product.

靶标

功能	The 26S protease is involved in the ATP-dependent degradation of ubiquitinated proteins. The regulatory (or ATPase) complex confers ATP dependency and substrate specificity to the 26S complex.
序列相似性	Belongs to the AAA ATPase family.
细胞定位	Cytoplasm. Nucleus.

图片



Western blot of proteasome 19S subunit S4 on HeLa cell lysate using ab3317.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors